

# Jazzy<sup>®</sup>

## ELITE HD

The Jazzy<sup>®</sup> Elite HD offers heavy-duty construction, a 204 kg (450 lbs.) weight capacity, and excellent performance in a stylish, highly manoeuvrable package.



**Pride**<sup>®</sup>  
Mobility Products  
Australia Pty. Ltd.

## Specifications

<b>Weight capacity</b>	204 kg (450 lbs.)*
<b>Maximum speed</b>	Up to 6.4 km/h (4 mph)
<b>Ground clearance</b>	7.6 cm (3") at bottom of base
<b>Turning radius</b>	61 cm (24") with front riggings
<b>Overall length</b>	90.2 cm (35.5") without front riggings
<b>Overall width</b>	61.6 cm (24.25")
<b>Drive wheels</b>	35.6 cm (14") solid
<b>Front wheels</b>	7.6 cm (3") solid anti-tips
<b>Rear wheels</b>	15.2 cm (6") solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	60 amp controller, Dynamic Shark
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 5 amp (standard)
<b>Per-charge range</b>	Up to 15.77 km (9.8 miles)
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	NF-22
<b>Battery weight</b>	17 kg (37 lbs.) each
<b>Base weight</b>	44 kg (98.5 lbs.)
<b>Standard Seat weight</b>	23 kg (51.5 lbs.) (high-back)
<b>Max. seat dimensions</b>	56 cm x 56 cm (22" x 22")

\* Maximum weight capacity is reduced to 400 lbs. (181 kg) on 46 cm (18") deep seat configurations.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.



**Pride**<sup>®</sup>  
 Mobility Products  
 Australia Pty. Ltd.  
 — LIVE YOUR BEST —

20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

(Tel) (03) 8770 9600 • (Fax) (03) 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)

## Features

- Two-motor, front-wheel drive
- 35.5 cm (14") knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- 5 amp, off-board charger

## Options

- High-back seat with solid seat pan
- Cup holder
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## Colour Choices



Jazzy  
Red



Jazzy  
Blue



Quality  
 Endorsed  
 Company  
 9001 QEC20029  
 SAI Global