



Narrowest complex rehab base on the market









Specifications

DRIVE WHEELS	35.6 cm (14") or 31.8 cm (12.5") pneumatic or flat free tyres
CASTER WHEELS	Front: 12.7 cm (5") flat free Rear: 15.2 cm (5.9") flat free
MAXIMUM SPEED ¹	Up to 10 km/h Up to 5.6 km/h w/iLevel activated
GROUND CLEARANCE	7.4 cm (2.9")
TURNING RADIUS ²	53.3 cm (20.9")
OVERALL LENGTH ²	82.5 cm (32.4") without front riggings
BASE WIDTH	55.2 cm (21.7")/35.6 cm (14") d/wheels 52.7 cm (20.7")/31.8 cm (12.5") d/wheels
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 30.5 cm (12") - 51 cm (20") Depth: 30.5 cm (12") - 51 cm (20")
TRU-Balance® 3 Power Tilt	Width: 30.5 cm (12") - 51 cm (20") Depth: 30.5 cm (12") - 51 cm (20")
TRU-Balance® 3 Tilt and iLevel®	Width: 30.5 cm (12") - 51 cm (20") Depth: 30.5 cm (12") - 51 cm (20")
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	40.6 cm (15.9") - 45.7 cm (17.9")
TRU-Balance® 3 Power Tilt 5	44.5 cm (17.5") - 49.5 cm (19.4")
TRU-Balance® 3 Tilt & iLevel® 5	44.5 cm (17.5") - 49.5 cm (19.4")
MANUAL RECLINE 6	Yes (TRU-Balance® Flex only)
BATTERY SIZE ³	55 AH (2)
BATTERY WEIGHT ⁷	55 AH: 18.1 kg each
AVAILABLE ELECTRONICS 8	120A Q-Logic 3 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6®
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT WITHOUT SEATING SYSTEM	57.3 kg

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 30.5 cm (12") of seat elevation at walking speed (up to 5.6 km/h)
- Fender lights and USB charger included

Options

- iLevel® option allows 30.5 cm (12") of seat elevation while at walking speed (5.6 km/h)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogy vent tray⁹
- Compact vent tray⁹
- 30.5 cm (12) x 30.5 cm (12") seat sizes with iLevel (SOF only)

Colours*



- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 0.3 cm with 45.7 cm x 45.7 cm seat, fixed towers, 136 kg weight capacity and solid tires (add 0.95 cm with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 1.4 kg based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.





