

PRIDE
NEXTGEN
MANUAL WHEELCHAIR SERIES

Model - Transit

(Below)

PMW-NGT



Model - Self Propelled

(Above)

PMW-NGSP

User

135 kg

Weight Capacity

Self-Propelled Model

15.5 kg

Unit Weight


Transit Model

13.7 kg

Unit Weight

The Next Gen Manual Wheelchair Series from Pride Mobility® is engineered to meet the needs of active users. These robust, lightweight self-propelled and transit models offer easy-to-use operation and outstanding performance.

Available Colours



Model - Self-Propelled

PMW-NGSP

Weight Capacity	135 kg (297.6 lbs.)
Overall Length	100 cm (39")
Overall Height	96 cm (38")
Overall Width	70 cm (27")
Seat Width	46 cm (18")
Seat Depth	40 cm (16")
Seat Back Height	40 cm (16")
Seat-to-Floor	52 cm (20")
Rear Wheel	60.96 cm (24") quick-release spoke wheel, solid tyres
Front Castor	20.32 cm x 3.17 cm (8" x 1.25" solid tyres)
Ground Clearance	2 cm - 7 cm (0.78" - 2.75") (adjustable anti-tip wheels)
Total Weight	15.5 kg (34.17 lbs.)

Model - Transit

PMW-NGT

Weight Capacity	135 kg (297.6 lbs.)
Overall Length	98 cm (38")
Overall Height	96 cm (38")
Overall Width	70 cm (27")
Seat Width	46 cm (18")
Seat Depth	40 cm (16")
Seat Back Height	40 cm (16")
Seat-to-Floor	52 cm (20")
Rear Wheel	30.48 cm (12"), solid tyres
Front Castor	20.32 cm x 3.17 cm (8" x 1.25" solid tyres)
Ground Clearance	3.5 cm (1.4")
Total Weight	13.7 kg (30.2 lbs.)

Features

- Quick-release 24" wheels
- Attendant brakes (with lock)
- Swing-away legrest with calf support
- Height-adjustable angled armrest
- Flip-down backrest
- Seat belt
- Adjustable anti-tip wheels



Features

- Lightweight, powder coated aluminium frame
- Attendant brakes (with lock)
- Swing-away legrest with calf support
- Height-adjustable angled armrest
- Flip-down backrest
- Seat belt

Warranty

- 12-month warranty (both models)