



Solutions For Life



POWERED BY  **coolcore®**

stealthproducts.com

SIMPLICITY G



Features:

- Contoured, high-density molded foam base
- Coccyx and sacral relief
- Moisture resistant and easy to clean Stealth-Tek incontinence cover.
- Sizes available - 16" - 24"W x 16" - 20"D

SPECTRUM FOAM SPP



Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Multi-layer, viscoelastic foam insert under at risk bony areas
- CoolCore® technology cover standard.
- Sizes available - 16" - 32"W x 16" - 28"D

SOLUTION SPP



Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, aqueous gel insert under at risk bony areas
- CoolCore® technology cover standard.
- Urethane liner for incontinence
- Sizes available - 10" - 24"W x 10" - 22"D

SPECTRUM GEL SPP



Features:

- Lateral and pre-ischial shelf contoured, high-density molded foam base
- Twin-cell, lightweight fluid insert under at risk bony areas
- CoolCore® technology cover standard.
- Sizes available - 14" - 24"W x 16" - 24"D

TRU-COMFORT 2 SPP



Convenient side pouch on both sides of the cushion with easy pull rubber loop zipper handles. Pouch is large enough to fit a larger format mobile phone.


Features:

- Deep contoured, high-density molded foam base
- CoolCore Technology cover standard.
- Easy pull rubber loop zipper handles
- Convenient side pouches on both sides of the cushion
- Sizes available - 16" - 22"W x 14" - 24"D

POWERED BY  **coolcore®**

What is CoolCore®

- Chemical Free material
- UPF/Flammability rated
- 4-Way Stretch nylon interlock with Elastane
- Patented combination of chemical-free blends that deliver 3 distinct functions:
Moisture Wicking to Moisture Transportation to Moisture Evaporation.
- **CoolCore® Technology is exclusive to Stealth Products.**



Year Round Temperature Control

CoolCore fabrics wick moisture away from the body in every direction keeping you dry and comfortable no matter what season it is.

Breakthrough Technology

Research driven, patented fabrics from Coolcore are chemical-free and yield remarkable benefits. They cool better, last much longer and offer protection from the sun's harmful rays (UPF).

Stealth Products Cushion Specifications

Cushion Name	Weight Capacity	Weight of Cushion	Cushion Height	Width Range (cm)	Depth Range (cm)
Simplicity G	158 kgs	997 g	5 cm	40.6 - 60.9	40.6 - 45.7
Solution SPP	158 kgs	2.4 kgs	7.6 cm	35.5 - 60.9	35.5 - 55.8
Spectrum Foam SPP	294 kgs	1.3 kgs	7.6 cm	40.6 - 81.2	40.6 - 71.1
Spectrum Gel SPP	158 kgs	1.5 kgs	7.6 cm	35.5 - 60.9	45.7 - 60.9
Tru Comfort 2 SPP	204 kgs	1.8 kgs	11.4 cm	40.6 - 55.8	40.6 - 60.9



Learn more at **StealthProducts.com**

QUANTUM®

www.pridemobility.com.au

www.quantumrehab.com



[/quantumrehab](https://www.facebook.com/quantumrehab)



[@quantum_rehab](https://www.instagram.com/quantum_rehab)