

Pride® Walkers

PMW 462

- 4-WHEELED
- ALUMINIUM
- LOOP BRAKE
- 6" WHEEL

FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 6" x 1" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard



TECHNICAL SPECIFICATIONS

| | |
|------------------|---------------------|
| Seat height: | 55 cm |
| Width: | 63.5 cm |
| Height: | 85-95 cm |
| Length: | 57 cm |
| Seat Width: | 35 cm |
| Seat Depth: | 31 cm |
| Weight: | 6.5 kg |
| Weight Capacity: | 114 kg |
| Colour: | Red/Black Patterned |

PMW 463

- 4-WHEELED
- ALUMINIUM
- LOOP BRAKE
- 8" WHEEL

FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 8" x 1.5" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard



TECHNICAL SPECIFICATIONS

| | |
|------------------|---------------------|
| Seat height: | 55 cm |
| Width: | 63.5 cm |
| Height: | 88-98 cm |
| Length: | 60 cm |
| Seat Width: | 35 cm |
| Seat Depth: | 31 cm |
| Weight: | 7 kg |
| Weight Capacity: | 114 kg |
| Colour: | Red/Black Patterned |

Pride®
Mobility Products
Australia Pty. Ltd.

PMW 6

- 4-WHEELED
- ALUMINIUM
- LOOP BRAKE
- 6" WHEEL

FEATURES

- Lightweight durable frame
- Detachable padded curved backrest
- Adjustable height
- Four 6" x 1" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard



TECHNICAL SPECIFICATIONS

| | |
|------------------|--------------|
| Seat height: | 52.8 cm |
| Width: | 61.5 cm |
| Height: | 77.7-93.2 cm |
| Length: | 62.4 cm |
| Seat Width: | 30.5 cm |
| Seat Depth: | 3.5 cm |
| Weight: | 6 kg |
| Weight Capacity: | 136 kg |
| Colour: | Red |

PMW 8

- 4-WHEELED
- ALUMINIUM
- LOOP BRAKE
- 8" WHEEL

FEATURES

- Lightweight durable frame
- Detachable padded curved backrest
- Adjustable height
- Four 8" x 1.5" wheels
- Bolt head-inner tube knob interface design, tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard



TECHNICAL SPECIFICATIONS

| | |
|------------------|--------------|
| Seat height: | 58.4 cm |
| Width: | 62.5 cm |
| Height: | 81.7-97.2 cm |
| Length: | 65.4 cm |
| Seat Width: | 30.5 cm |
| Seat Depth: | 3.5 cm |
| Weight: | 6 kg |
| Weight Capacity: | 136 kg |
| Colour: | Red |



Quality
Endorsed
Company
BIO 9991 05209923
244 134-04

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

Pride[®]
Mobility Products
Australia Pty. Ltd.

20-24 Apollo Drive, Hallam Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

Tel 03 8770 9600 • Fax 03 9703 2720

www.pridemobility.com.au