

MANUAL WHEELCHAIR RANGE

Pride manual wheelchairs are designed to meet the needs of active users with feature-rich, stylish models. Pride manual wheelchairs are extremely lightweight, highly durable and provide enhanced fit and function.



Pictured is the L400. Not available in red.

L400

Overall Width	66 cm (26")
Overall Height	97 cm (38")
Overall Length	84 cm (33")
Folded Width	30 cm (11")
Total Weight	13 kg (29 lbs.)
Weight Capacity	113 kg (250 lbs.)
Seat Dimensions	46 cm x 40 cm (18" x 16")
Seat-to-floor Height	47 cm - 50 cm (19" - 20")
Seat Back Height	38 cm - 48 cm (15" - 19")
Colour	Silver

FEATURES

- Adjustable tension seat upholstery increases comfort and positioning
- Flip-up, desk length armrests improve accessibility
- Multiple-position rear wheel and caster mounts optimise centre-of-gravity and seat-to-floor height adjustments
- Quick-release, removable rear wheels
- Full-profile 61 cm x 3.5 cm rear wheels for enhanced outdoor capabilities
- Long-handle wheel locks ensure comfort and security
- Through-bolted, industry-compatible leg rests with aluminium foot plates for durability

PMW900



FEATURES

- Easy-to-use folding system ideal for car travel or compact storage, including a flip down back rest
- Flip-up armrests and adjustable and detachable leg rests
- Large self-propelling, 24" pneumatic or solid tyres with quick release wheels and sturdy 8" solid aluminium casters
- Comes standard with a calf strap and seat belt
- Available in your choice of solid or pneumatic tyre models. Pneumatic tyre model comes standard with tyre pump.
- Available in 'silver liquid' colour

SPECIFICATIONS

Overall Width	78 cm (30.7")
Overall Height	93 cm (36.6")
Overall Length	103 cm (40.5")
Folded Width	26.5 cm (10.5")
Total Weight	21 kg (46.3 lbs)
Weight Capacity	110 kg (242 lbs)
Weight w/out Leg rest & Wheels	12.2 kg (27 lbs)
Seat Dimensions	45 cm x 40 cm (18" x 15.7")
Seat-to-floor Height	50 cm (19.7")
Seat Back Height	46 cm (18.1")

PMW901



FEATURES

- Ideal economy transit chair with 12" solid rear wheels and solid 8" front aluminium casters
- Convenient detachable, swing away leg rests improves manoeuvring as well as portability
- Fixed backrest and fixed armrests for added stability and support
- Features a rear back pocket for users to carry belongings
- Available in 'silver liquid' colour

SPECIFICATIONS

Overall Width	59.5 cm (23.4")
Overall Height	91 cm (35.8")
Overall Length	96 cm (37.8")
Folded Width	25 cm (9.8")
Total Weight	15 kg (33 lbs)
Weight Capacity	113 kg (249 lbs)
Weight w/out Leg rest & Wheels	Unable to remove wheels
Seat Dimensions	46 cm x 40 cm (18" x 15.7")
Seat-to-floor Height	50 cm (19.7")
Seat Back Height	41 cm (16.1")

PMW902



FEATURES

- Large self-propelling, 24" pneumatic or solid tyres with quick release wheels and sturdy 8" solid aluminium casters
- Designed to be lightweight and easily transported
- Useful flip-down back rest for compact storage
- Anti-tip wheels for added safety and security
- Comes standard with a calf strap and seat belt
- Flip-up armrests and adjustable and detachable leg rests
- Available in your choice of solid or pneumatic tyre models. Pneumatic tyre model comes standard with tyre pump.
- Available in: *Metallic Blue* colour

SPECIFICATIONS

Overall Width	65 cm (25.6")
Overall Height	95 cm (37.4")
Overall Length	103 cm (40.5")
Folded Width	30 cm (11.8")
Total Weight	13.8 kg (30.4 lbs)
Weight Capacity	110 kg (242 lbs)
Weight w/out Leg rest & Wheels	9.5 kg (21 lbs)
Seat Dimensions	46 cm x 41 cm (18" x 16")
Seat-to-floor Height	49 cm (19.3")
Seat Back Height	47 cm (18.5")

PMW903



FEATURES

- Ideal economy transit chair with 12" solid rear wheels
- Sturdy aluminium construction
- Features height adjustable anti-tip wheels for added safety
- Comes standard with a calf strap and seat belt
- Easy-to-use folding system ideal for car travel or compact storage, including a flip down back rest
- Flip-up armrests improve both accessibility and manoeuvring in and out of the chair
- Available in: *Metallic Blue* colour

SPECIFICATIONS

Overall Width	65 cm (25.6")
Overall Height	95 cm (37.4")
Overall Length	96 cm (37.8")
Folded Width	30 cm (11.8")
Total Weight	12.5 kg (27.5 lbs)
Weight Capacity	110 kg (242 lbs)
Weight w/out Leg rest & Wheels	Unable to remove wheels
Seat Dimensions	46 cm x 41 cm (18" x 16")
Seat-to-floor Height	49 cm (19.3")
Seat Back Height	47 cm (18.5")

PMW904



FEATURES

- Ideal lightweight economy transit chair with 8" solid wheels
- Fixed arm rests
- Useful flip-down back rest for compact storage or travel
- Comes standard with a calf strap and seat belt
- Features adjustable and detachable leg rests
- Available in: *Red Metallic Liquid* colour

SPECIFICATIONS

Overall Width	57.8 cm (22.7")
Overall Height	96 cm (37.8")
Overall Length	94 cm (37")
Folded Width	25.4 cm (10")
Total Weight	11 kg (24.2 lbs)
Weight Capacity	110 kg (242 lbs)
Weight w/out Leg rest & Wheels	Unable to remove wheels
Seat Dimensions	48 cm x 40 cm (18.9" x 15.7")
Seat-to-floor Height	48 cm (18.9")
Seat Back Height	47 cm (18.5")

PMW905



FEATURES

- Ideal lightweight economy aluminium transit chair with solid 12" rear wheels
- Fixed arm rests
- Comes standard with a calf strap and seat belt
- Features easy-to-use loop brakes
- Useful flip-down back rest for compact storage or travel
- Available in: *Red Metallic Liquid* colour

SPECIFICATIONS

Overall Width	57.8 cm (22.7")
Overall Height	96 cm (37.8")
Overall Length	94 cm (37")
Folded Width	25.4 cm (10")
Total Weight	11 kg (24.2 lbs)
Weight Capacity	110 kg (242 lbs)
Weight w/out Leg rest & Wheels	Unable to remove wheels
Seat Dimensions	48 cm x 40 cm (18.9" x 15.7")
Seat-to-floor Height	48 cm (18.9")
Seat Back Height	47 cm (18.5")