

# POWERASSIST



The PowerAssist features an easy-to-use squeeze trigger mechanism along with a speed dial, enabling users to customise speed settings.



## POWER UP

The Pride Power Assist is designed to aid caregivers in pushing a manual wheelchair. The PowerAssist attaches directly on to manual wheelchairs easily.

**Pride**<sup>®</sup>  
Mobility Products  
Australia Pty. Ltd.

— LIVE YOUR BEST —

20-24 Apollo Drive, Hallam VIC 3803  
PH: (03) 8770 9600 | FAX: (03) 9703 2720

[www.pridemobility.com.au](http://www.pridemobility.com.au)  
ACN 088 609 661 ABN 66 088 609 6611



# POWERASSIST

## SPECIFICATIONS

Maximum Speed<sup>1</sup>: Up to 5 km/h (3.1 mph)

Range per Charge<sup>1</sup>: Up to 10 km/h (6.2 mph)

Maximum Gradient: 6°

Weight Capacity: 120 kg (264 lbs)

Motor: 12V 100W 4100RPM

Total Weight: 14 kg (30.8 lbs)

Battery: 12V 17ah Sealed Lead Acid  
Battery Weight: 6.05 kg (13.3 lbs)

Charger: 3A 12V output

Fits Wheelchair Widths: 16" - 20"

Dimensions<sup>2</sup>:  
Height: 285mm (11.2")  
Length: 280mm (11")  
Width: 220mm (8.6")



*The Power Assist Battery Pack*



*The Power Assist Unit*

<sup>1</sup> Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

<sup>2</sup> Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

## FEATURES

- The Pride Power Assist is designed to aid caregivers in pushing a manual wheelchair in otherwise strenuous situations such as up steep hills, inclines and loose gravel.
- The Power Assist is designed to work in conjunction with most wheelchairs with seat widths of 16"-20".
- Weight capacity of 120 kg/264 lbs (wheelchair and user combined).
- The Power Assist has a maximum speed of up to 5 km/h (3.1 mph).

**Pride**<sup>®</sup>  
Mobility Products  
Australia Pty. Ltd.  
— LIVE YOUR BEST! —

20-24 Apollo Drive, Hallam VIC 3803  
PH: (03) 8770 9600 | FAX: (03) 9703 2720

www.pridemobility.com.au  
ACN 088 609 661 ABN 66 088 609 6611