

SPORT RIDER

**3 WHEEL
SCOOTER**



Pride[®]
Mobility Products
Australia Pty. Ltd.
— LIVE YOUR BEST™ —

Tel 03 8770 9600 • Fax 03 9703 2720
www.pridemobility.com.au

SPECIFICATIONS

Maximum Speed ¹	Up to 15 kph
Ground Clearance ²	13cm (5.1")
Front Tyre	35.6cm (14")
Rear Tyre	35.6cm (14")
Cruise Range	45 - 50 km
Turning Radius	157.3cm (62")
Incline Capabilities	10°
Battery	70AH MK Gel Batteries x 2
Charger Capacity	8A Off-board
Motor Output	3.5 HP

DIMENSIONS


Length ²	165cm (65")
Width ²	78cm (30.7")
Height	112cm (44")
Wheel Base	110cm (43.6")
Seat Width	50.4cm (20")
Seat Depth	50cm (19.8")
Back Height	73cm (28.9")
Seat Height from Deck	47.1cm (18.7")
Seat Height from Ground	55cm (21.8")

WEIGHT

Base Section	98kg (215.6lb)
Battery (x 2)	68kg (149.6lb)
Seat	22.5kg (49.5lb)
Total Weight with Battery ³	188.5kg (417.7lb)
Weight Capacity	180kg (400lb)

AVAILABLE COLOURS

Black



Warranty 12 months

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³Total weight is measured with maximum 100amp batteries.

NOTE: All specifications subject to change without notice.

Pride[®]
 Mobility Products
 Australia Pty. Ltd.

20-24 Apollo Drive, Hallam Victoria 3803
 ACN 088 609 661 ABN 66 088 609 661
 Tel 03 8770 9600 • Fax 03 9703 2720
 www.pridemobility.com.au



FEATURES

- Full digital displaying speedometer/odometer
- Full suspension system
- High back, adjustable seat
- Powerful front/rear lighting kit
- Large rear storage pod



*Sport rider
 with rear pod*

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

PRIDE[®]
 PERFORMANCE MOBILITY
 — LIVE YOUR BEST[®] —

