



30.5 cm x 30.5 cm
(12" x 12")
seat sizes now
available with
iLevel®

Q6 edge[®] 2.0

With features like standard 9.6 km/h (6 mph) motors and advanced responsiveness, the Q6 Edge[®] 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel[®] technology.

quantumrehab.com

AU-06ED6E2-021516

20-24 Apollo Drive
Hallam, Victoria 3803
www.pridemobility.com.au
ACN 088 609 661 ABN 66 088 609 661



Tel: 03 8770 9600
Fax: 03 9703 2720

The Rehab Power Chair
Solutions Company

Specifications

DRIVE WHEELS	Standard: 35.56 cm (14") Solid Optional: 35.56 cm (14") Pneumatic
CASTOR WHEELS	Front: 15.24 cm (6") Solid (standard) 20 cm (8") Solid (optional) Rear: 15.24 cm (6") Solid (standard) 20 cm (8") Solid (optional)
	N/A
MAXIMUM SPEED ⁴	Up to 9.6 km/h (6 mph)
GROUND CLEARANCE	6.9 cm (2.7") (battery tray)
TURNING RADIUS ⁵	52.1 cm (20.5")
OVERALL LENGTH	90.2 cm (35.5") without foot riggings
BASE WIDTH	62.2 cm (24.4")
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 25 cm - 56 cm (10-22") Depth: 25 cm - 56 cm (10-22")
TRU-Balance® 3 Power Tilt	Width: 30.5 cm - 56 cm (12-22") Depth: 30.5 cm - 56 cm (12-22")
TRU-Balance® 3 Lift & Tilt	Width: 41 cm - 56 cm (16-22") Depth: 35.5 cm - 56 cm (14-22")
SEAT-TO-FLOOR HEIGHTS	
Synergy®/TRU-Balance® 3 Static Seat ⁸	42 cm - 47 cm (16.5-18.5")
TRU-Balance® 3 Power Tilt ⁸	44 cm - 49.5 cm (17.5-19.5")
TRU-Balance® 3 Lift & Tilt ⁸	44 cm - 49.5 cm (17.5-19.5")
MANUAL TILT/RECLINE ¹⁰	Yes
BATTERY SIZE ⁶	NF-22, 55 AH
BATTERY WEIGHT ¹⁴	17.24 kg (38 lbs.)
AVAILABLE ELECTRONICS ¹²	100A Q-Logic 2 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6®
WEIGHT CAPACITY ⁷	136 kg (300 lbs.) ⁹
BASE WEIGHT	63 kg (138.6 lbs.)


- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- .3175 cm (1/8") with 46 cm x 46 cm (18" x 18") seat, fixed towers, 136 kg (300 lbs.) weight capacity and solid tyres (add .95 cm (3/8") with pneumatic tyres). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 113 kg (250 lbs.) weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 34 kg (75 lbs.) and maximum dimensions are 35.5 cm x 35.5 cm (14" W x 14") D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 1.4 kg (3 lbs.) based on manufacturer.
- 15) Call for availability.

Features

- Standard 9.6 km/h (6 mph) motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- iLevel® option available with up to 25.4 cm (10") of seat elevation (in 16 seconds) at walking speed (up to 5.6 km/h (3.5 mph))
- iLevel is retrofittable on the Q6 Edge® 2.0

Options

- iLevel® option allows 25.4 cm (10") of seat elevation in 16 seconds while at walking speed (5.6 km/h (3.5 mph))
- TRU-Balance® 3 Seating
- TRU-Comfort 2 Seating System
- TRU-Balance® 3 25.4 cm (10") Power Adjustable Seat with iLevel
- Q-Logic 2 EX Drive Controls
- Trilogy vent tray¹⁵
- Compact vent tray¹⁵
- LED full lighting package
- 10 cm (4") wide flat-free knobby extended traction tyres¹⁵
- 30.5 cm x 30.5 cm (12" x 12") seat sizes with iLevel (SOF only)
- 20 cm (8") solid castors
- New LED Fender Lights



NEW

NEW LED FENDER LIGHTS

Colours

SEAT BACKS for units without TRU-Balance® 3 are only available in Back in Black.



(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • 20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

www.pridemobility.com.au • www.quantumrehab.com • www.ilevel.rehab

