

XLR USB CHARGER



Shown with cell phone holder



Plugs conveniently into
the charger ports of
mobility scooters and
power chairs

Features and Benefits

- Charge your cell phone while on the go
- Can be used while power chair or scooter is in motion
- Charging detection—switches off automatically
- Standard with lanyard
- Short circuit, thermal and under-voltage protection
- XLR USB Charger will not fully discharge the battery

Compatible with: Most iPhone®, Android and mobile devices

*iPhone® is a registered trademark of Apple, Inc.

5V DC, 1 A charger that is compatible with 24V and 36V systems

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Jazzy
POWER CHAIRS

GO GO®
TRAVEL MOBILITY

pridemobility.com.au



20-24 Apollo Drive Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661
(Tel) 03 8770 9600 • (Fax) 03 9703 2720

