



QUANTUM[®]

quantumrehab.com





T H E 4 F R O N T O F A N E W E R A

4FRONT

The 4Front™ is a quiet, more responsive front-wheel drive power chair that features automotive-grade suspension with unprecedented comfort and ride quality. The 4Front's CASE (Caster Angle Sensor Encoder) technology and SRS (Smooth Ride Suspension) provide stability for smooth obstacle transitions and unmatched driving performance. The 4Front™ features standard speeds of up to 10 km/h in the lowered position, and speeds up to 5.1 km/h while elevated at 10" with our optional safe seat elevation. Optional LED fender lights allow you to see and be seen.

QUANTUM® 4FRONT shown in Pretty in Pink with 10" Safe Seat Elevation and Q-Logic 3 Drive Controls



5.1 KM/H
**SAFE SEAT
ELEVATION**
(OPTIONAL)

FRONT AND REAR
LED FENDER LIGHTS

CASE (CASTER ANGLE
SENSOR ENCODER)
TECHNOLOGY

SRS
(SMOOTH RIDE SUSPENSION)







iLevel Power Chairs feature a revolutionary new power chair design for a better experience. iLevel improves daily functionality and quality of life by allowing operation of the power chair with the seat elevated while at walking speed (up to 5.6 km/h).

12" of power adjustable seat height while at walking speed (10" on Q6 Edge® HD)

iLevel's patented Extra Stability Technology™ enhances safety for transfers, reach, and numerous other activities of daily living

12 exciting colour options

Retrofittable on the Q6 Edge® 2.0

136 kg weight capacity on Q6 Edge® 2.0 and Q6 Edge® Z

181 kg weight capacity on Q6 Edge® HD

Extra Stability Technology™ for a smooth and stable ride

Front and rear LED fender lights are optional on the Q6 Edge® 2.0, Q6 Edge® Z and Q6 Edge® HD with iLevel®

12" LIFT
AT
5.6 KM/H

Maximum lift of 10" on Q6 Edge HD



Q6 SERIES

power bases



Q6 Edge 2.0

With features like standard 10 km/h motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

FEATURES

- User weight capacity of 136 kgs.
- Mid-Wheel 6® allows six wheels on the ground for stability
- Standard 10 km/h motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- iLevel® option available with up to 12" of power adjustable seat height at walking speed (5.6 km/h)
- Available with Q-Logic 3 Advanced Drive Control System
- 12 exciting colour options
- Available with WC-19 crash-test approved occupied transit kit option on tilt, tilt & recline and iLevel®
- Available with Accu-Trac Advanced Tracking Technology
- Optional switch-operated front and rear LED fender lights
- Optional 8" casters (no additional charge)

Q6 Edge Z

The Q6 Edge® Z utilises powerful 4-pole motors for maximum torque and Mid-Wheel 6® design for aggressive outdoor performance and precision indoor manoeuvrability. The Q6 Edge® Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability, and is available with our iLevel® technology.

FEATURES

- User weight capacity of 136 kgs
- Mid-Wheel 6® allows six wheels on the ground for stability
- 10 km/h motor package option available
- iLevel® option available with up to 12" of power adjustable seat height at walking speed (5.6 km/h)
- 12 exciting colour options
- Optional switch-operated front and rear LED fender lights
- Available with Accu-Trac Advanced Tracking Technology





Quantum 610

The compact Quantum® 610 is the ultimate choice for users who demand the superb stability of patented Mid-Wheel 6® Design and tight-quarter manoeuvrability. Highly adaptable, the Quantum 610 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.

FEATURES

- User weight capacity of 136 kgs
- Mid-Wheel 6® allows six wheels on the ground for stability
- Patented ATX Suspension (Active-Trac® with extra stability) for enhanced performance over more varied terrain
- Compatible with TRU-Balance® Power Positioning
- Maximum speed of 6.4 km/h
- 10" drive wheels

Jazzy 623

With features like standard 4-pole motors, Mid-Wheel 6® drive design and ATX Suspension, the Jazzy® 623 is engineered to meet the needs of the most active user. The Jazzy 623 is available in several different seating options, including highback seat and full rehab configuration. The Dynamic Shark controller is easy to operate and offers two power seating functions.

FEATURES

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with REAC lift and tilt system as well as TRU-Balance® Power Positioning



Jazzy 623 in tilt position

- Low impact OMNI-Casters (nylon, spherical-shaped castors) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI- Castors and semi-independent rear castor beam for enhanced performance over more varied terrain



Jazzy 623 with black ADI-back seat

REAR-WHEEL DRIVE

power bases



Quantum R44 (Rival)

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space manoeuvrability. The Rival accepts a complete range of rehab seating and electronics options.

FEATURES

- 7 exciting colour options
- Compact size for outstanding manoeuvring in tight spaces
- Maximum speed up to 10 km/h with high-speed motor package
- High-performance, high efficiency motor package
- Accepts a wide range of seating and positioning systems
- Two suspension packages
- Compatible with TRU-Balance® 3 Power Positioning Systems
- Sliding battery tray for easy access

PEDIATRIC

power bases

Quantum 4Front

The 4Front™ features standard speeds of up to 10 km/h in the lowered position, and speeds up to 5.1 km/h while elevated at 10" with our optional safe seat elevation. Optional LED fender lights allow kids to see and be seen.

FEATURES

- iLevel® now available for seat sizes as small as 12" x 12"
- CASE (Caster Angle Sensor Encoder) technology
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Superior, automotive-grade SRS
- Available with optional 10" of poweradjustable safe seat elevation at speeds of up to 5.1 km/h
- Q-Logic 3 EX Advanced Drive Controls
- Optional LED front and rear fender lights
- Bold aesthetics with 13 colour options



Q6 Edge 2.0

With features like standard 10 km/h motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the needs of the most active child.

The Q6 Edge® 2.0 is available with our iLevel® technology and features bold aesthetics with 12 colour options to fit any child's personality. It accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

FEATURES

- iLevel® now available for seat sizes as small as 12" x 12"
- Standard 10 km/h motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics with 12 colour options
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 5.6 km/h)
- iLevel® is retrofittable on the Q6 Edge® 2.0
- Available with Q-Logic 3 Advanced Drive Control System
- Available with WC-19 crash-test approved occupied transit kit option on title, tilt & recline and iLevel
- Optional front and rear LED fender lights
- Available with Accu-Trac Advanced Tracking
- Optional 8" casters (no additional cost)
- Three suspension packages available



BARIATRIC

power bases



Q6 Edge HD

With features like standard 4-pole motors, Mid-Wheel 6® Drive Design and ATX Suspension, the Q6 Edge® HD is engineered to meet the performance needs of the most active user. The Q6 Edge HD accepts a complete range of seating and electronics options to deliver exceptional rehab capability.

FEATURES

- 204 kg weight capacity (181 kgs with iLevel® option)
- iLevel® option available with up to 10" of power adjustable seat height at walking speed (5.6 km/h)
- 12 exciting colour options
- Mid-Wheel 6® for maximum stability
- Compatible with TRU-Balance® 3 HD Power Positioning
- ATX Suspension (Active-Trac® with extra stability)
incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- Optional switch-operated front and rear LED fender lights


Quantum 1450

The Quantum® 1450 front-wheel drive bariatric power base features innovative design that provides exceptional outdoor performance while delivering outstanding tight space manoeuvrability. The 1450 accepts a complete range of rehab seating and electronics options.

FEATURES

- 272 kg weight capacity
- 272 kg heavy duty tilt capacity
- Max. speed up to 7.2 km/h
- High-performance, high-efficiency motor package
- Compatible with TRU-Balance® Power Positioning and a wide range of other seating and positioning systems





TRU-Balance® 3 Power Positioning Systems

- Make fine adjustments while client is in the chair
- Order your seat width independent of your back width at no additional charge
- Internally or externally rotate cantilever armrests for optimal positioning
- Easily retrofit power tilt and power recline
- iLevel® is retrofittable on the Q6 Edge® 2.0
- Easy back height and depth adjustment
- Easy center of gravity adjustment
- A full range of premium accessories available
- Optional attendant control



TRU-BALANCE® 3

power positioning systems



Power Recline

- 85° - 175° of recline
- 136 kg weight capacity
- Does not increase seat-to-floor height
- Standard shear reduction and optional cantilever armrests
- Available in up to 12 back shroud colours
- Available with TRU-Balance® 3 Power Tilt, TRU-Balance® 3 Power Adjustable Seat Lift and iLevel.



Power Tilt

- 50° of tilt
- Lower standard seat-to-floor height
- Unmatched adjustment and durability
- Maintains center of gravity shift
- Available in up to 12 back shroud colours



Power Adjustable Seat Lift with iLevel®

- iLevel® technology enables the unit to drive up to 5.6 km/h while fully elevated
- Offers 12" of power adjustable seat lift (10" maximum lift on Q6 Edge® HD)
- Available in up to 12 back shroud colours
- Available on the Q6 Edge® 2.0, Q6 Edge® Z, and Q6 Edge® HD



Power Articulating Foot Platform

- Broad articulation range for enhanced adjustment
- 204 kg weight capacity
- Multiple calf pads, calf panel and footplate options
- Unmatched adjustment, durability and strength



Heavy Duty Power Tilt & Recline

- 204 kg weight capacity
- 50° of tilt
- 85° - 175° of recline
- Available in up to 12 back shroud colours

Please Note:

Some options may not be available on all power bases

SEATING

systems



Q6 Edge® 2.0 shown
with Synergy® Seating
System

Synergy® Seating System

The Quantum® Synergy® Seating System offers a clean, sleek appearance and unparalleled ease of adjustment. The seating system's unique design allows for height, width, depth and back angle adjustments to be made quickly and easily in the field. The Synergy Seat interfaces with most types of seat and back cushions, and is designed to be the most user-friendly system available today.

FEATURES

- Easily adjustable and versatile design that can grow with the client
- Easy customization
- Several width/depth adjustment packages ranging from 10" to 24"
- Accepts a wide variety of after-market seats and backs
- Heavy Duty option provides a weight capacity of 272 kgs and seating up to 32" x 28" (model specific)

ADI Q-Back

The ADI Q-Back from Stealth Products, Inc. integrates seamlessly with the TRU-Balance® 3 Seating Systems and is offered in three configurations:

- Active Contour ADI Q-Back (2")
- Deep Contour ADI Q-Back (4")
- Extra Deep Contour ADI Q-Back (6")

Available sizes for tilt and static seating:

Heights: 19"-22"
Widths: 16", 18", 20"

Available sizes for power recline seating:

Heights: 22"-25"
Widths: 16", 18", 20"

Note: Additional back sizes are available with Special Order.




Active Contour
ADI Q-Back



Deep Contour
ADI Q-Back



Extra Deep Contour
ADI Q-Back



Our ADI catalogue provides you with the very latest product information, right at your fingertips.



Stealth Products Cushions Colour Guide

stealthproducts.com



Synergy® Simplicity

Tremendous comfort
and moderate pressure
management



FEATURES

- High-density polyurethane foam construction
- Exclusive wrap-and-lift design for pressure management
- Four-way stretch, fluid resistant upper cover
- Coccyx and sacral relief



Synergy® Spectrum®

Customized comfort and superb
positioning and pressure
management



FEATURES

- Extra deep well
- Choice of twin-cell lightweight fluid or viscoelastic foam insert
- Viscoelastic top layer
- 2-cover system to prevent shearing with interior liner
- 294 kg weight capacity with viscoelastic foam insert

POWERED BY  coolcore®



Synergy® Solution®

Three-component design for
pressure management, positioning
and long-lasting comfort



FEATURES

- Exclusive wrap-and-lift design for pressure management
- 2-cover system to prevent shearing with interior liner
- Twin-cell gel insert
- Viscoelastic top layer



TRU-Comfort 2 Back and Seat Cushion



FEATURES

- Adjustable positioning and pressure management
- Compatible with both TRU-Balance® 3 Power Positioning Systems and the Static Seat Frame System
- Available with four-way stretch, recovery, fluid-resistant, breathable cover
- Additional 1" of viscoelastic foam for maximum comfort and pressure management

ELECTRONIC

drive controls

Q-LOGIC 3

Advanced Drive Control System

Featuring



Q-Logic 3 EX
Joystick Controller



Q-Logic 3 EX
Enhanced Display

Q-Logic 3 Drive Control System

The Q-Logic 3 Advanced Drive Control System from Quantum Rehab® is like no other power wheelchair electronics system. Designed from clinician and consumer feedback, Q-Logic 3 is user-friendly and highly customisable. This adaptable, expandable control system can increase and maximise independence.



Home Screen



Drive Screen

FEATURES

- Wireless Bluetooth programming through Windows devices
- Bluetooth access to computer, tablet, and phone functions
- iAccess programmable push button and toggle options for customised seating access
- Simple firmware updates
- Capable of utilising multiple input devices
- Preset drive settings
- Preset colour background options
- Clinic Mode:
 - Quick setup of multiple input devices during client evaluations
 - Automatically assigned drive modes to input devices

Q-Logic 3 EX Enhanced Display

The optional Q-Logic 3 EX Enhanced Display serves as a separate screen used with specialty drive controls to display system data. With the Enhanced Display's built-in environmental controls, Q-Logic 3 can operate a wide range of appliances that use infrared (IR) remote controls, such as televisions, and serves as a Bluetooth® computer mouse emulator.



Optional Modules



**Q-Logic 3 SAJ/
Attendant Control
Joystick**



Q-Logic 3 iAccess
iAccess programmable
push button and toggle
options (not pictured)
for customised seating
access



**Q-Logic 3 Specialty Controls
Interface Module**

ELECTRONIC

drive controls

Non-Expandable



Dynamic Shark



PG VR2 (4-key)



PG VR2 (6-key)



PG VR2 (9-key)

Expandable (Q-Logic 3)



Q-Logic 3 EX
Joystick



Q-Logic 3 EX
Enhanced
Display



Q-Logic 3
EX Stand Alone
Joystick / Attendant
Joystick



Q-Logic 3
Specialty
Controls
Interface
Module

Expandable (Q-Logic 2)



Q-Logic 2 EX
Joystick



Q-Logic 2 EX
Enhanced
Display




Q-Logic 2
EX Stand Alone
Joystick /
Attendant
Joystick



Q-Logic 2
Specialty
Controls
Interface
Module



Stealth® i-Drive®
Tri-Array Series,
3-Switch



Stealth® Ultra
3-Switch Head
Array

Proportional Drive Controls



Mo-Vis™
All-Round
Joystick



Mo-Vis™
All-Round Light
Joystick



Mo-Vis™
Multi Joystick



Mo-Vis™
Micro Joystick



Stealth®
Precision Mini
Proportional
Joystick



Stealth®
Mushroom
Joystick

Single Switch



Stealth®
Egg Switch
(Multiple colours
available)



Mo-Vis™
Twister Basic
Button



Mo-Vis™
Twister Satellite
Button

Sip-N-Puff



Stealth®
Sip-N-Puff


Misc. Options



Joystick
Handles

TRU-BALANCE® 3

accessories



Accessory Bar



Transfer Bar
(Pair-Quick Release)



Mobile Phone
Holder



Glove Box



Cup Holder
(Generic)



Personal Item
Hook



Backpack Holder (Pair)



iPad Holder



Oxygen Holder



Cup Holder
























Reflector Kit



Push Handles (Pair)

COLOUR

options


		Edge 2.0	4Front	J623	Edge Z	Q610	Edge HD	R44	Q1450
Deep Purple		•	•		•		•		
Midnight Blue		•	•		•		•		
Viper Blue				•		•		•	
Electric Blue		•	•		•		•		
True Blue								•	
Emerald Green								•	
* Night Train			•						
Onyx Black						•		•	
Back in Black		•	•		•		•		
Grey Street		•							
Silver Bullet		•							
White		•							
Lime Light		•	•		•		•		
Yellow Blaze		•							
Orion Yellow								•	
Sunburst Orange						•			
Orange Crush		•							
Candy Apple Red				•		•		•	•
Red Corvette		•	•		•		•		
Pearl Pink						•		•	
Pretty in Pink		•							

*Available on 4Front only

Please note optional colours are available upon special request.


SERVICE

and support




quantumrehab.com

The ultimate resource for everything you need to know about Quantum® products.



Clinical Resources

A special section of **quantumrehab.com** for clinicians, featuring clinical references and more.



Follow us...

Stay up to date on the latest Quantum® events and product releases through social media.



stealthproducts.com

Stealth Products® launched in 1999, in response to a need for versatile, high-quality head and neck positioning systems, as well as switch site location hardware. Their products include: head and body support/positioning, joystick mounts, trays, seating, and pediatric mobility products.

Quantum Rehab® Professional App

The Quantum Rehab® Professional App is available for Apple and Android smart phones and tablets. Look for the app in the Apple App Store or Google Play links using your smart device to download the app directly.



Sales and Technical Support

For more informaton, contact your local Authorised Quantum Provider



3/66 Pilkington St, Garbutt QLD 4814

wheelzmobility@hotmail.com.au

0422 106 511

Stamp/Quantum Provider Details



OPTIONAL:
10" LIFT AT
5.1 KM/H

	Q6 EDGE [®] 2.0	AVAILABLE WITH iLEVEL [®]	4FRONT [™]	AVAILABLE WITH 10" LIFT	J623	Q6 EDGE [®] Z	AVAILABLE WITH iLEVEL [®]
DRIVE CONFIGURATION	Mid-wheel Drive		Front-wheel Drive		Mid-wheel Drive		Mid-wheel Drive
DRIVE WHEELS	14" Flat Free/ 14" Pneumatic		14" Flat Free/ 14" Pneumatic		14" Flat Free/ 14" Pneumatic		14" Flat Free/ 14" Pneumatic
CASTER WHEELS	Front 6" Solid (Optional 8") Rear 6" Solid (Optional 8")		Rear: 9" Flat Free		Front 6" Solid Rear 6" Solid		Front 6" Solid Rear 6" Solid
ANTI-TIP WHEELS	N/A		Front: 4"		N/A		N/A
MAXIMUM SPEED ⁴	Up to 10 km/h		Up to 10 km/h		Up to 8.6 km/h		Up to 10 km/h
GROUND CLEARANCE	2.7" (battery tray)		3" (motor pivot)		3" (motor pivot)		3" (frame base)
TURNING RADIUS ⁵	6" casters: 20.5" 8" casters: 23.25"		24.5"		20"		19.75"
OVERALL LENGTH ¹¹	6" casters: 35.5" 8" castors: 40.48"		45.7"		35.5"		35.6"
BASE WIDTH	24.4"		24.2"		24"		26.25"
Seating Sizes:	TRU-Balance [®] 3		TRU-Balance [®] 3		TRU-Balance [®] 2		TRU-Balance [®] 3
SYNERGY [®]	W: 10-22" D: 10-22"		N/A		W: 14-22" D: 12-22"		W: 10-24" D: 10-24"
STATIC	W: 10-22" D: 10-22"		W: 12-22" D: 12-22"		W: 10-22" D: 10-22"		W: 12-24" D: 12-24"
POWER TILT	W: 12-22" D: 12-22"		W: 12-22" D: 12-22"		W: 14-22" D: 12-22"		W: 12-24" D: 14-24"
LIFT & TILT	W: 12-22" D: 12-22"		W: 12-22" D: 12-22"		W: 12-22" D: 12-22"		W: 12-24" D: 12-24"
Seat-to-Floor Heights: ¹⁶	TRU-Balance [®] 3				TRU-Balance [®] 2		TRU-Balance [®] 3
SYNERGY [®] / STATIC ⁸	16.5-18.5"		15.25-17.25" <small>NOT AVAILABLE WITH SYNERGY[®]</small>		15.9-17.9"		16.5-18.5"
POWER TILT ⁸	17.5-19.5"		17.5-19.5"		17.3-19.3"		17.5-19.5"
LIFT & TILT ⁸	17.5-19.5"		17.75-19.75"		17.5-19.5"		18-20"
MANUAL RECLINE	Yes		Yes		Yes		Yes
BATTERY SIZE ⁶	NF-22 (55Ah)		Group 24: 23.3 kgs (55Ah) Group 34: 17.4 kgs (60Ah)		NF-22 (55Ah)		Group 24 (75Ah)
BATTERY WEIGHT ¹⁴	17.2 kgs		23.3 kgs		17 kgs		23.5 kgs (each)
AVAILABLE ELECTRONICS ¹²	120A Q-Logic 3 EX		120A Q-Logic 3 EX		75A Dynamic Shark		100A Q-Logic 2 EX
BATTERY CHARGER	8A Off-Board		8A Off-Board		8A Off-Board		8A Off-Board
MOTOR PACKAGES	2-Pole Mid-Wheel 6 [°]		2-Pole front-wheel drive		4-Pole mid-wheel		4-Pole mid-wheel
WEIGHT CAPACITY ⁷	136 kgs		136 kgs		136 kgs		136 kgs
BASE WEIGHT	62.8 kgs (6" casters) 67.1 kgs (8" casters)		63.5 kgs		62 kgs		65 kgs

1) Standard

2) Optional

3) Requires HD Suspension

4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

5) Without front riggings

6) Battery capacity may change when adding power positioning systems.

7) Weight capacity is dictated by the selection of seating.

8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 136 kgs weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in their lowest position; does not include elevating system unless specified.

9) 113 kgs weight capacity with TRU-Balance[®] Systems

10) A minimum seating size of 20" x 20" is required with iLevel[®] on the Q6 Edge[®] HD; maximum seating size is 24" x 25"

11) Without foot riggings

12) All Quantum units can be ordered with specialty controls

13) iLevel[®] not available with 12.8 km/h motor package

14) Battery weight may vary +/- 1.3 kgs based on manufacturer.

15) Weight capacity is reduced to 181 kgs with iLevel[®] power adjustable seat height

16) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.



QUANTUM® 610	Q6 EDGE® HD		QUANTUM® R44® (RIVAL)		QUANTUM® 1450
Mid-wheel Drive	Mid-wheel Drive		Rear-wheel Drive		Front-wheel Drive
10" Flat Free/ 10" Pneumatic	14" Flat Free/ 14" Pneumatic		14" Flat Free/ 14" Pneumatic		14" Flat Free
Front 5" Solid Rear 6" Solid	Front 6" Solid Rear 6" Solid		Front: 9" Solid		9" Solid
N/A	N/A		N/A		4" Solid
Up to 6.4 km/h	Up to 7.5 km/h		Up to 9.6 km/h		Up to 7.2 km/h
2.5" (frame)	3" (frame centre)		3.25" (frame base)		2" (anti-tip wheels)
22"	19.75"		35.75"		26.5"
34.74"	35.6"		36"		38.25"
23.19"	26.25"		25"		29.25"
TRU-Balance® 2	TRU-Balance® 2 HD	TRU-Balance® 3 HD	TRU-Balance® 2 HD	TRU-Balance® 3 HD	TRU-Balance® 2 HD Tilt
W: 10-20" D: 10-20"	W: 20-28" D: 20-28"	W: 20-32" D: 20-28"	W: 10-22" D: 10-20"	W: 12-22" D: 12-22"	W: 20-32" D: 20-28"
W: 10-20" D: 10-20"	W: 20-28" D: 20-28"	W: 20-32" D: 20-28"	W: 10-22" D: 10-20"	W: 12-22" D: 12-22"	W: 20-32" D: 20-28"
W: 14-20" D: 14-20"	W: 20-28" D: 20-28"	W: 20-32" D: 20-25"	W: 14-22" D: 14-20"	W: 12-22" D: 12-22"	W: 20-32" D: 20-28"
W: 16-20" D: 16-20"	N/A	W: 20-24" D: 20-25"	N/A	W: 12-22" D: 12-22"	N/A
16.5-17.5"	16.5-18.5"	17.5-19.5"	16.5-18"	16.5-18.1"	16.5-18.5"
16.875"	18.2-20.9"	18.5-20.5"	16.5-18"	18-19.6"	18.125-20.125"
17.75"	N/A	17.75-19.75"	N/A	18.5-20.1"	N/A
Yes	No		Yes		No
U-1 (31Ah)	Group 24 (75Ah)		Group 24 (75Ah)		Group 24 (75Ah)
10.6 kgs	23.5 kgs		23.5 kgs		23.5 kgs
100A Q-Logic 2 EX	100A Q-Logic 2 EX		90A VR2 6-key 100A Q-Logic 2 EX		90A VR2 6-key 100A Q-Logic 2 EX
8A Off-Board	8A Off-Board		8A Off-Board		8A Off-Board
2-Pole Glide 6.4 km/h ¹	4-Pole 7.5 km/h ¹		4-Pole 10 km/h ¹		4-Pole 7.2 km/h
113 kgs ⁹	204.1 kgs ¹⁵		136 kgs		272 kgs
45.4 kgs	65 kgs		58.5 kgs		88.4 kgs

The information contained herein is correct at the time of publication.
We reserve the right to alter specifications without prior notice.



(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • 20-24 Apollo Drive Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

www.pridemobility.com.au • www.quantumrehab.com



/quantumrehab



@quantum_rehab